Original

THE SYNTHETIC & ART SILK MILLS' RESEARCH ASSOCIATION

(Linked to the Ministry of Textiles, Govt. of India) Sasmira Marg, Worli, Mumbai - 400 030.

Tel.: 2493 5351 / 52, 2493 2047, 2493 8753 / 54 (EPBX) • Fax: 91-22-2493 0225 • Grams: "SASMIRA" BOMBAY-DADAR 400 030. E-mail: ed@sasmira.org • testing@sasmira.org • Web: www.sasmira.org

SASMIRA TEST LABORATORIES

TEST REPORT

Test Report No.:	SASMIDA/TD 20/10-20		
Name & Address of customer Sample forwarding letter No.	SASMIRA/TR-30/19-20	Date:	
	Clique Medical Solutions Pvt. Ltd., 2B, 1806, Dreams Complex, Behind Dreams Mall, Bhandup (W), Mumbai – 400 078		
	Letter dated 04.04.19		
Particulars of sample given by the customer	Sample ID Mark		
Date of receipt	1. Blue Colour Jacket (Thermaissance)		ce)
Date of testing	04.04.19		
	04.04.19 to 05.04.19		

Name of the Test : Autoclaving cycle

Additional Test Information:

a) Autoclaving Temperature: 121°C

b) Autoclaving Time: 20 mins

Test method	D
Autoclave cycle	No change observed after one Autoclaving cycle at 121°C for 20 m

Note: This test is not under NABL scope of accreditation

Authorized Signatory

Bhakti Ghag Microbiology Lab

-- End of the report --

Note: 1. The result relate only to the samples submitted by customer for testing

The test report shall not be reproduced except in full, without written approval

of the laboratory and shall not be used for publicity or litigation. 3. The sample/s tested and reported is/are not drawn by Laboratory.

Catalyze your business in association with sasmira